

Before use, clean tube with alcohol wipe and/or gauze

NOT PROVIDONE—IODINE

1. Place the white neckband behind the patient's neck ensuring tabs are facing outwards close to the corners of the mouth.

2. For fitting of patients neck size the neckband can be adjusted by removing the double sided Velcro tab, folding excess material underneath and refixing the velcro tab in the correct position.

3. Make sure the tube is cleaned with an appropriate wipe and when dry, remove adhesive hook tab from the blue securing strap. Peel off paper backing and adhere around the top of the endotracheal tube at the lip/teeth line.

4. Centralise the blue securing strap on top of the tube, covering the adhesive hook tab using the wider central section of the blue strap.

5. Wrap both sides of the blue securing straps around the tube for a complete rotation, crossing at the top. Ensure the adhesive hook tab is completely covered.



6. Twist the ends of the blue securing strap 180 degrees around each other to reduce/eliminate any slipping on the tube. (Do not tie a knot)

7. Fasten the end sections of the blue securing straps onto the white velcro tabs attached to the neckband ensuring sufficient tension to stabilise the tube.

8. The tautness can be adjusted by lifting the double sided velcro tabs while keeping the blue strap ends attached and relocating to tighten the strap and ensure security of the tube.

9. It is recommended these tabs are checked for secureness after being repositioned. Once secured, the blue securing strap should only be removed from the velcro tabs when they are being replaced. New blue strap can then be attached.

10. Replace blue securing strap when wet or soiled, with spare strap provided. Replace white neckband when required.

Upon removal of the tube holder, please ensure that all three parts of the holder are removed safely, accounted for and disposed of in the correct manner.

