



All Silicone 2 way Foley Catheter

Description

All silicone indwelling Foley 2 way catheter, traditional balloon construction, hand finished to meet exact surgical and patient requirements. BS 1616.

Sizes

GMS Code	Description	Reference	GMS Price
93317	Male 5-10ml 14fg	CAT84M-14/SI	
93441	Male 30ml 18fg	CAT85M-18/SI	
93441	Male 30ml 16fg	CAT85M-16/SI	
93317	Male 5-10ml 20fg	CAT84M-20/SI	
93317	Male 5-10ml 22fg	CAT84M-22/SI	
93317	Male 5-10ml 12fg	CAT84M-12/SI	
93317	Male 5-10ml 16fg	CAT84M-16/SI	
93317	Male 5-10ml 18fg	CAT84M-18/SI	
93317	Male 5-10ml 24fg	CAT84M-24/SI	
93838	Female 30ml 20fg	CAT87F-20/SI	
93838	Female 30ml 18fg	CAT87F-18/SI	
93838	Female 30ml 16fg	CAT87F-16/SI	
93843	Female 5-10ml 22fg	CAT86F-22/SI	
93843	Female 5-10ml 20fg	CAT86F-20/SI	
93843	Female 5-10ml 18fg	CAT86F-18/SI	
93843	Female 5-10ml 16fg	CAT86F-16/SI	
93843	Female 5-10ml 14fg	CAT86F-14/SI	
93843	Female 5-10ml 12fg	CAT86F-12/SI	
93891	Paediatric 3-5ml 10fg	CAT83P-10/SI	
93891	Paediatric 3-5ml 8fg	CAT83P-8/SI	



Indication

For use when either short or medium term catheterisation is required.

Advantages

- Manufactured from high quality Silicone.
- Established design (Foley type).
- Silicone elastomer equates to low reactivity by body tissues, minimised urethral irritations and reduced catheterisation complications, such as trigonitis, urethritis and meatitis.
- Non-wetting surface facilitates drainage, is resistant to encrustation by urinary salts and is non-adherent, providing for greater patient comfort.
- Fewer catheter changes mean that infection caused by repeated removal is markedly reduced.
- Translucent with blue radio-opaque line.
- All our product range is latex and phthalate free.