




**INSIGHT**  
Medical Products Ltd

PRODUCT CATALOGUE 2015



[www.insightmedical.net](http://www.insightmedical.net)  
+44 (0)1666 500055




# AN INSIGHT INTO INSIGHT

Insight Medical Products Ltd was established by it's current directors in May 2000. From our base in the Cotswolds, Insight has become a highly respected supplier of good value, high quality innovative products to the healthcare market.

We are pleased to work with healthcare professionals to design and manufacture innovative products that will help budgets and improve patients well-being while receiving treatment.

Insight now has it's own in-house manufacturing facility with a fully accredited cleanroom facility for manufacturing, assembly and packaging of our products. All of our devices comply to the Medical Device Directive and our regulatory systems meet with the requirements of EN ISO 13485:2012 and Directive 93/42/EEC under the requirements supervision of SGS.

We pride ourselves on our reputation for Customer Service and are committed to ensuring a first class service to all of our customers, in the UK and overseas.



# CONTENTS


## TRACHEOSTOMY CARE

- Tracheostomy Tube Holder
- Tracheostomy Tube Cleaning Swabs
- Advazorb Plus Dressing
- Advazorb Lite Dressing
- Advazorb Fixation Dressing



## TUBE FIXATION

- Insil Silicone Tape
- Endotracheal Tube Holder
- Transducer Holder
- Tube Tidy Clips



## LARYNGECTOMY PROTECTORS

- Sofnex Laryngectomy Protectors
- Sofshield Laryngectomy Protectors
- Cascade Protectors




## LARYNGECTOMY CARE

- Advabond Silicone Glue
- Purifoam Laryngectomy Filters
- Actibalm Manuka Honey Skincare
- Voice Prosthesis Brushes




## AIRWAY MANAGEMENT

Tracheal Tube Introducers  
Guedel Airways



## THEATRE CONSUMABLES

Post Operative Eye Shields  
Insight Bite Guard



## MISCELLANEOUS

Curve-Hook  
Intermittent Catheters



## VETERINARY

Cat Bladder Irrigation Catheters  
Dog Bladder Irrigation Catheters



# TRACHEOSTOMY CARE

## TRACHEOSTOMY TUBE HOLDER

CE

The Insight Tracheostomy Tube Holder uses the highest quality materials to provide a safe and secure fixation for the tracheostomy tube.

The neck band is manufactured from an ultra soft brushed cotton/foam material that prevents soreness and discomfort to the patient.

A simple to use two-piece design allows quick fitting to the majority of neck sizes, while the integral 'cough-strip' enables reflex movement to the patient if needed.

- Adult and Paediatric sizes available
- Available on FP10 Prescription
- 2-Piece design allows secure fitting for most sizes
- Cough reflex section offers more patient comfort
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
TH/100	Tracheostomy Tube Holder	Adult	289-4814	FTH000	10
TH/150	Tracheostomy Tube Holder	Paediatric	289-4806	FTH092	10

## TRACHEOSTOMY TUBE CLEANING SWABS




The Insight Cleaning Swabs are perfect for cleaning all types of PVC and silicone tracheostomy or laryngectomy tubes.

Available in two sizes, our Tube Cleaning Swabs are secure and suitable for popular sizes of tubes up to 6mm (small) or 9mm (large).

The soft foam heads of the swabs avoid damaging or scuffing the inside of the tube, whilst the white plastic handle is flexible for ease of use.

- Use with warm water to remove mucus residue in the tube.
- Available on FP10 Prescription
- Two available sizes suitable for popular sizes of tubes
- Latex-Free.
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
TCS/900	Tube Cleaning Swab	Large	324-8689	FDH420	100
TCS/950	Tube Cleaning Swab	Small	324-8671	FDH421	100

SIZE	HEAD WIDTH	HEAD LENGTH	HANDLE LENGTH	MAX TUBE
LARGE	7/16"	1"	6"	9mm
SMALL	1/4"	7/8"	6"	6mm

# TRACHEOSTOMY CARE

## ADVAZORB PLUS (T) TRACHEOSTOMY DRESSING

Advazorb Plus (T) offers a new standard in hydrophilic foam tracheostomy dressings with incredible absorbency ability.

Manufactured in a lightweight, soft and conforming low-adherent foam with high fluid absorbance and retention characteristics.

Newly improved with a film backing added to ensure improved strength and absorbance for high/medium exuding wounds.

Supplied sterile and individually wrapped, ready for use.

Available in one convenient size to suit most intubation procedures, and to offer comfort to the patient.


- Cushioning effect on wound
- Low adherence – Reduced trauma on removal
- Highly absorbent 5mm foam
- Unbeatable Comfort
- Rapid uptake of exudates
- Can be easily cut to size
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



CE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3900	Advazorb Plus Dressing	80mm x 80mm	-	FDH419	20





## ADVAZORB (T) LITE TRACHEOSTOMY DRESSING

Advazorb (T) Lite is the definitive lightweight hydrophilic foam tracheostomy dressing.


Advazorb (T) offers remarkable standards of conformability, softness, lightness and air permeability.

Perfect for the treatment of light/medium exuding wounds, as moisture is absorbed and retained in the hydrophilic cellular structure. No loose fibres.

Ultra efficient performance.

Supplied sterile and individually wrapped, ready for use. Available in one convenient size to suit most intubation procedures, and to offer comfort to the patient.

- Cushioning effect on wound
- Low adherence – Reduced trauma on removal
- Highly absorbent 2.5mm foam
- Available on FP10 Prescription
- Unbeatable comfort
- Rapid uptake of exudates
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/4212	Advazorb Lite Dressing	80mm x 80mm	374-0230	ELY392	20

# TRACHEOSTOMY CARE

## ADVAZORB FIXATION TRACHEOSTOMY DRESSING

CE

Advazorb Fixation is a post decannulation tracheostomy dressing fitted with a central location button that prevents unwanted escape of air through the wound site with gentle pressure applied during occasions such as talking, coughing and sneezing.


Advazorb Fixation is a hydrophilic foam dressing with a unique perforated soft silicone wound contact layer and surrounding adhesive. It is an atraumatic foam dressing which does not adhere to moist wounds.

The silicone helps maintain a moist wound healing environment and makes it ideal as a non-adherent wound contact dressing. Soft silicone dressings do not damage newly formed granulating tissue and do not lose their adhesion after initial application.

The silicone border prevents the need for additional dressing or fixation. Silicone adhesive combines security with pain-free dressing changes. Waterproof backing enables the patient to shower without the need to replace or cover the dressing.

Supplied sterile and individually wrapped, ready for use.

- Location button enables accurate identification of the tracheostomy site
- Supplied sterile, for single patient use
- Maintains a moist environment and manages exudate
- Allows atraumatic dressing change
- Provides a barrier for bacteria
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/4341	Advazorb Fixation	12cm x 10cm	-	FDH605	10

## INSIL SILICONE TAPE



This versatile tape is backed with a silicone adhesive. When tape is removed it is extremely gentle on the skin, however it can be re-applied without losing adhesive qualities.

Insil multi-use tape can be used to fix most medical devices such as tubes, IV cannulae, drains, probes, eye-shields and dressings.

Used in theatres for securing closure of eyelids during operations, and allowing painless removal.

Insil is available in 2cm and 4cm widths with fine perforations every 5cm to allow simple fixing without using scissors. It is supplied in pocket-sized dispenser boxes for ease of use.

- Available on FP10 Prescription
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
INSIL/200	Insil Silicone Tape	20mm x 3m	342-4132	-	12
INSIL/400	Insil Silicone Tape	40mm x 1.5m	342-4140	-	12

# TUBE FIXATION

## ENDOTRACHEAL TUBE HOLDER



The Insight Endotracheal Tube Holder Set has been commended by professionals on the ease of use, security and patient comfort.

Because of the soft materials, the products has been proven to greatly reduce the incidents of sores to the patients mouth, while maintaining a secure fitting to the ultra soft neckpiece.

Each set is supplied with 3 holders to enable the user to replace soiled items, therefore improving value over alternative methods.

- Fully adjustable to fit most neck sizes
- Laminated hook/loop tabs offer total security of fitting
- Ultra soft materials prevent abrasions to the patient
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	NHS CODE	PACK QTY
TFP/200	Endotracheal Tube Holder Set	Up to 20" neck	FTH005	10
TFP/201	Endotracheal Tube Holder Set Extra Long	Up to 24" neck	FTH025	10
TFP/222	Endotracheal Tube Holder Neck Band	Up to 20" neck	-	10

## HALLAMSHIRE TRANSDUCER HOLDER




A simple and effective method for ensuring that the transducer remains in the correct position either within the department or while the patient is being moved and replaces the need for adhesive tape application.

The holder ensures stable positioning of the transducer for extended periods of time whilst providing patient comfort.

Made from soft material and secured to the arm with fully adjustable hook and loop fastener to fit most patients, eliminating the difficulties in taping the transducer to the patient or repositioning the transfusion stand due to patient movement, also saving nursing time.

- Designed in partnership with the ICU at Royal Hallamshire Hospital in Sheffield
- For use in the department or while patient is being transferred
- No need for any use of adhesive tape on the patient
- Quality soft foam/cotton material for patient comfort
- Easy to apply
- Latex-Free
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
TFP/400	Transducer Holder	Standard	-	-	10

# TUBE FIXATION

## TUBE TIDY CLIPS

CE

The simplest product on the market for keeping tubes and lines arranged and sorted during patient monitoring.

The product can clip to blankets, sheets, pillows and opens to allow tubing to run through.

Simply lay tubes required over hook/loop and seal to hold tubes in place.

Organises cords and drapes, useful in all specialist units and wards.

- Can retain drapes and linens without the use of pins or fixers, thus eliminating the chance of personal injury
- User is able to write on the product for identity markings
- Simply open the hook/loop material, insert tubing and close for use
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
TFP/800	Tube Tidy Clips	White	-	-	10
TFP/850.B	Tube Tidy Clips	Blue	-	-	20
TFP/850.R	Tube Tidy Clips	Red	-	-	20
TFP/850.Y	Tube Tidy Clips	Yellow	-	-	20

# LARYNGECTOMY PROTECTORS

## SOFNEX LARYNGECTOMY PROTECTORS



Sofnex Protectors can be washed a minimum of 10 times, offering a significant cost saving over alternative protectors that still use foam inserts, and reduces visits to the doctor for new prescriptions.

Sofnex is available in 3 colours, White, Blue or Beige. The blue and beige protectors have an antibacterial coating to help prevent infections.

Sofnex can be used with cravats or other decorative covers, although the appearance of Sofnex allows the user to wear it without a cover.

Sofnex is available in one size only, and is supplied with two neckstraps in every box of ten protectors.

- Available on FP10 Prescription
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3856	Sofnex Protector	White	319-0402	-	10
CR/3859	Sofnex Protector	Blue	319-0410	-	10
CR/3860	Sofnex Protector	Beige	319-0428	-	10

# LARYNGECTOMY PROTECTORS

## SOFSHIELD LARYNGECTOMY PROTECTORS



Sofshield Laryngectomy Protectors can be washed a minimum of 10 times, offering a significant cost saving over alternative protectors that still use foam inserts.

Available in 4 colours, white, blue, beige and black. The blue and beige protectors have an antibacterial coating to help prevent infections occurring, while also eliminating the need for costly additional cosmetic covers.

Sofshield uses a unique toggle system that eliminates the need to tie cords.

- Available on FP10 Prescription
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
AS/3669	Sofshield -White	Large	291-9231	FTH006	10
AS/3764	Sofshield -Blue	Large	305-5373	FTH008	10
AS/3765	Sofshield -Beig	Large	305-5381	FTH007	10
AS/3698	Sofshield -White	Small	291-9249	FTH009	10
AS/3777	Sofshield -Blue	Small	305-5399	FTH100	10
AS/3778	Sofshield -Beige	Small	305-5407	FTH010	10
AS/4051	Sofshield -Black	Medium	377-1722	-	10



## CASCADE SHOWER PROTECTOR



Cascade is a soft waterproof PVC protector that provides a barrier to prevent shower water from entering a stoma thereby improving the patient's confidence to shower routinely.

The neckband is soft and elastic to form a comfortable yet firm seal around the neck to prevent water seeping below the protector. The backing is a knitted polyamide material that will absorb moisture from surrounding wet skin and provide a comfortable warm feel, unlike other products that may initially feel cold to the skin.

From a safety perspective, the polyamide backing also prevents adhesion to the skin or stoma button when bib is wet, as adhesion from unlined protectors can reduce the ability for a patient to breath properly.

For added comfort and breathability the Cascade Protector can be worn over a Sofshield Protector or an alternative filter, however, this is not essential due to the backing material. The protector is an aesthetically pleasing white in colour and supplied in pairs to allow for drying between showers.



- Available on FP10 Prescription
- Shower proof
- Lay flat contour
- Washable
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
AS/3835	Cascade Protector	White	314-2601	FTH028	10 (5 packs of 2)

# LARYNGECTOMY CARE

## ADVABOND SILICONE GLUE

CE

Advabond Silicone Glue offers the user a high grade, skin friendly silicone adhesive for attaching stoma appliances to the skin to reinforce fixation of an adhesive base plate to skin around the stoma.

Advabond ensures that the appliance remains adhered to the skin even in the presence of moisture or perspiration due to it's high viscosity.

Advabond has long lasting non-gumming and non-peeling characteristics.

- Apply direct to appliance or onto skin
- Available on FP10 Prescription
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
AS/3896	Advabond Silicone Glue	40ml	328-6413	-	1

## PURIFOAM LARYNGECTOMY FILTERS



Purifoam Filters use a high quality medical grade foam that protects the stoma from dirt and dust entering the airway, while warming and humidifying the incoming air.

Purifoam is adhered to the skin using a skin friendly adhesive strip that has rounded edges for added comfort.

Manufactured from beige coloured foam to give a natural aesthetic look. Purifoam has excellent moisture retaining qualities, ensuring high levels of absorbency to give exceptional humidification to the airway for breathing reassurance.

- Available on FP10 Prescription
- Change daily or more frequent if required
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3853	Purifoam	6.5 x 6.5cm	319-0386	FTH029	30

# LARYNGECTOMY CARE

## ACTIBALM MANUKA HONEY SKINCARE

CE

Actibalm is a topical ointment made up from only two ingredients, white pharmaceutical grade petroleum jelly and 100% Activon Manuka Honey.

The petroleum jelly seals the skin preventing bacteria from entering the affected area and keeps skin moisture from evaporating whilst the honey's antibacterial properties can help heal and soothe.

Using Activon Manuka Honey, blended with petroleum jelly, not only protects the stoma, but heals skin simultaneously ensuring a healing process to tender or cracked skin and stomas.

Activon Honey is a medical grade Manuka Honey which has a reliable level of antibacterial potency not found in other honey products.

Activon honey is filtered, sterilised and controlled to ensure effective and safe application in wound management.

- Anti-bacterial , effectively kills harmful bacteria
- Maintains moist wound healing environment
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3955	Actibalm	10g	-	-	1


## VOICE PROSTHESIS BRUSHES



Voice Prosthesis Brushes are manufactured with soft nylon bristles that prevent damaging the silicone inner surface of the prosthesis.

Brushes are tapered to provide effective cleaning of mucus and secretion build-up in the inner lumen of the voice prosthesis.

- Available on FP10 Prescription
- Silicone safety sleeve to prevent over-insertion
- Flexible plastic handle for ease of use
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
PB/50	Voice Prosthesis Brush	Standard	229-4395	-	10


## TRACHEAL TUBE INTRODUCERS



Insight Tracheal Tube Introducers are designed to help insert tracheal Tubes during difficult intubation.

All introducers are single use to reduce chances of cross-infection and are produced in straight or angled format to suit all procedures.

- Sheen surface gives less tracheal tube resistance
- Graduated markings to enhance entry and exit management
- Barium tip to assist X-Ray control
- Available with straight or angled tip
- Full range of lengths and sizes available
- Saves costs in cleaning and product deterioration
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	SIZE	LENGTH	NHS CODE	PACK QTY
SB/05/ST	Sterile Bougie –Straight	5ch	470mm	FTH031	10
SB/10/600/S	Sterile Bougie –Straight	10ch	600mm	FTH032	10
SB/10/800/S	Sterile Bougie –Straight	10ch	800mm	FTH033	10
SB/15/600/S	Sterile Bougie –Straight	15ch	600mm	FTH034	10
SB/15/800/S	Sterile Bougie –Straight	15ch	800mm	FTH035	10
SB/10/600/A	Sterile Bougie –Angled	10ch	600mm	FTH036	10
SB/10/800/A	Sterile Bougie –Angled	10ch	800mm	FTH037	10
SB/15/600/A	Sterile Bougie –Angled	15ch	600mm	FTH040	10
SB/15/800/A	Sterile Bougie –Angled	15ch	800mm	FTH041	10
SBV/14/750/OX	Angled, Vented with Oxy. Con	14ch	750mm	FTH042	10

## GUEDEL AIRWAYS




Insight high quality Guedel Airways are available in eight various sizes, with different colour coded rigid bite blocks to ensure ease of identification for use.

Our Guedel Airways are CE marked, latex free and provided as clinically clean bulk packed.

The Airways smooth lumen dimension allows highly competent patient aspiration.

- Manufactured from high quality materials
- Colour coded for ease of identification
- Full range of sizes
- Rigid Bite Block
- Smooth lumen and external surface for easy insertion
- Rounded atraumatic edges
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE




IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
OA/50/000	Clear	000	-	-	50
OA/50/00	Blue	00	-	-	50
OA/50/0	Black	0	-	-	50
OA/50/1	White	1	-	-	50
OA/50/2	Green	2	-	-	50
OA/50/3	Yellow	3	-	-	50
OA/50/4	Red	4	-	-	50
OA/50/5	Orange	5	-	-	50

## POST OPERATIVE EYE SHIELDS



Eye Shields are used after surgery to protect the eye and to hold dressings in place while holes in the shield provide ventilation.

Supplied in bulk for trade suppliers of Theatre Procedure Packs. We also supply individually sterile packed for hospitals and can use customers own artwork if required.



- Universal fitting available for dual fitting
- Available in Clear, Blue Tint, or Opaque White
- Left and Right fitting available
- UK Manufactured
- Latex-Free
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	COLOUR	PACK QTY
ES/500/C	Bulk Packed Eye Shield	Clear	2000
ES/500/W	Bulk Packed Eye Shield	White	2000
ES/500/B	Bulk Packed Eye Shield	Blue Tint	2000
ES/100/C	Sterile Packed Eye Shield	Clear	100
ES/100/W	Sterile Packed Eye Shield	White	100
ES/100/B	Sterile Packed Eye Shield	Blue Tint	100
ES/500/WR	Bulk Packed Eye Shield – Right	White	2000
ES/500/WL	Bulk Packed Eye Shield – Left	White	2000
ES/100/WR	Sterile Packed Eye Shield – Right	White	100
ES/100/WL	Sterile Packed Eye Shield – Left	White	100




## INSIGHT BITE GUARD



Insight Bite Guards are a protective shield used to protect teeth against accidental damage from contact with rigid instruments during surgery.

Enhanced design to provide safer retention when fitted. Manufactured from medical grade TPE for pliability and patient comfort.

Coloured bright yellow to allow operator visibility.



- Shaped for optimum fit to patient.
- Radio-opaque, autoclavable on the 134° cycle up to 3.5 minutes
- Can be supplied sterile or unsterile in easy peel blister pack
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	COLOUR	SIZE	NHS CODE	PACK QTY
BG/100	Unsterile Bite Guard	Yellow	Adult	FTH124	20
BG/200	Sterile Bite Guard	Yellow	Adult	FTH125	20

## CURVED AMNIOTIC HOOKS




Curve-Hook Anatomically Shaped Amniotic Membrane Perforator.

The Curve-Hook is anatomically shaped to make the hook more versatile and easier to use. The Curve-Hook has been designed to make the membrane easily accessible.

The design ensures optimum efficiency and the prevention of injury to the patient when used.

The shaped rounded tip and protected hook simply engages and ruptures the membrane, saving the need to prod or puncture, therefore reducing any chance of injury to the foetus.

- Individually supplied sterile in easy to open packs.
- Latex-Free
- UK Manufactured
- CE Marked



LATEX  
FREE



IMP CODE	DESCRIPTION	NHS CODE	PACK QTY
KH/100	Curve-Hook	-	100

## HANDHELD INTERMITTENT CATHETERS




The Handheld Intermittent Catheter is fully adjustable for comfort and ease of use for self catheterisation.

Movable handle to ascertain the ideal fitting position.

The catheter has a nelaton tip and made from PVC to minimise allergic reactions

- Individually supplied sterile in easy to peel pouches
- Latex-Free
- CE Marked



LATEX  
FREE

IMP CODE	HANDLE SIZE	CATHETER SIZE	PACK QTY
INS/100/8	115mm Length	180mm x 8ch	100
INS/100/10	115mm Length	180mm x 10ch	100
INS/100/12	115mm Length	180mm x 12ch	100
INS/100/14	115mm Length	180mm x 14ch	100

## CAT BLADDER IRRIGATION CATHETERS

Manufactured as an OEM product for the veterinary trade only, we do not sell direct to vets.

Single use disposable catheter with smooth closed tip and two opposed side eyes, all with adjustable suture collar and welded luer connector.

Available in two materials, polypropylene (without stylet) or polyethylene (complete with introducing stylet), for longer term treatment.

- Supplied packaged to customers requirements
- Blister packed, sterile, with customers own artwork, boxes of 12 with labelling
- Bulk packed, unsterile for customers own packaging
- Latex-Free
- UK Manufactured



LATEX  
FREE



IMP CODE	CATHETER SIZE	MATERIAL	STYLET	LUER
CC/03/130P	3FG 1.0mm x130mm Length	Polypropylene	No	Pink
CCS/03/130P	3FG 1.0mm x130mm Length	Polyethylene	Yes	Pink
CC/03/160P	3FG 1.0mm x160mm Length	Polypropylene	No	Pink
CCS/03/160P	3FG 1.0mm x160mm Length	Polyethylene	Yes	Pink
CC/04/130W	4FG 1.0mm x130mm Length	Polypropylene	No	White
CCS/04/130W	4FG 1.0mm x130mm Length	Polyethylene	Yes	White
CC/04/160W	4FG 1.0mm x160mm Length	Polypropylene	No	White
CCS/04/160W	4FG 1.0mm x160mm Length	Polyethylene	Yes	White


## DOG BLADDER IRRIGATION CATHETERS

Manufactured as an OEM product for the veterinary trade only, we do not sell direct to vets.

Single use disposable catheter with smooth closed tip and two opposed side eyes and welded luer connector.

Available in polypropylene (without stylet)

- Supplied packaged to customers requirements
- Blister packed, sterile, with customers own artwork, boxes of 12 with labelling
- Bulk packed, unsterile for customers own packaging
- Latex-Free
- UK Manufactured



LATEX  
FREE



IMP CODE	CATHETER SIZE	MATERIAL	STYLET	LUER
DC/04/500/W	4FG 1.3mm x 500mm	Polypropylene	No	White
DC/06/500/Y	6FG 2.0mm x 500mm	Polypropylene	No	Yellow
DC/08/500/O	8FG 2.6mm x 500m	Polypropylene	No	Orange
DC/10/500/R	10FG 3.3mm x 500mm	Polypropylene	No	Red







Units 3, Priory Industrial Estate  
Tetbury  
Gloucestershire  
United Kingdom GL8 8HZ  
Tel: +44 (0) 1666 500055  
Fax: +44 (0) 1666 500115  
Email: [info@insightmedical.net](mailto:info@insightmedical.net)



[Facebook.com/insightmedical](https://www.facebook.com/insightmedical)  
[@InsightMP](https://twitter.com/InsightMP)



[www.insightmedical.net](http://www.insightmedical.net)

Registered in England No: 3949101