



STAINLESS STEEL TROLLEYS GENERAL FEATURES

Totally made of 18/10 stainless steel, bright polish finishing. Heavy duty structure. Easy cleaning.

With 2 handles, shaped trays.

Frame :

18/10 stainless steel tube diam. 25 mm, castors diam. 75 mm.

2 MODELS :

Knockdown construction or argon welded.

KNOCKDOWN CONSTRUCTION 2 OR 3 TRAYS TROLLEYS

Shaped trays. Height from floor to trays :
2 trays : 300 and 820 mm. - 3 trays : 200, 510 and 820 mm
Bumpers diam. 95 mm.

2 TRAYS

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
441	600 X 430 mm	445	600 X 430 mm
741	750 X 500 mm	745	750 X 500 mm
941	900 X 600 mm	945	900 X 600 mm

3 TRAYS

WITHOUT GUARD RAILS		WITH GUARD RAILS	
REF.	SIZE	REF.	SIZE
446	600 X 430 mm	447	600 X 430 mm
746	750 X 500 mm	747	750 X 500 mm
946	900 X 600 mm	947	900 X 600 mm

3 TRAYS DRESSING TROLLEYS

EQUIPED WITH GUARD RAILS AND STAINLESS STEEL WASTE BIN 12 L

3 removable and height adjustable rails.

Removable waste bin holder, left and right directionable. Overall size 1080x570 mm

REF.	SIZE	CASTORS
752	750 X 500 mm	100 mm
752S	Argon welded version	

2 TRAYS TROLLEYS ARGON WELDED

WITHOUT GUARD RAILS WITH GUARD RAILS

REF.	SIZE	REF.	SIZE
441S	600 X 430 mm	445S	600 X 430 mm
741S	750 X 500 mm	745S	750 X 500 mm
941S	900 X 600 mm	945S	900 X 600 mm

3 TRAYS TROLLEYS ARGON WELDED

WITHOUT GUARD RAILS WITH GUARD RAILS

REF.	SIZE	REF.	SIZE
446S	600 X 430 mm	447S	600 X 430 mm
746S	750 X 500 mm	747S	750 X 500 mm
946S	900 X 600 mm	947S	900 X 600 mm

GUARD RAILS

REF.	SIZE
441 E6	for length 600 mm
741 E6	for length 750 mm
941 E6	for length 900 mm
441 E5	for width 400 mm
741 E5	for width 500 mm
941 E5	for width 600 mm
P27	clip for guard rail

TROLLEYS

TRAYS TROLLEYS ARGON WELDED TURNED DOWN RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø100mm. Delivered without accessories.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
451T	600 X 400 mm	456T	600 X 400 mm
751T	700 X 480 mm	756T	700 X 480 mm



TRAYS TROLLEYS ARGON WELDED RAISED RIM SHELVES

Made of 18/10 stainless steel, argon welded
Castor Ø100mm. Delivered without accessories.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
451R	600 X 400 mm	456R	600 X 400 mm
751R	700 X 480 mm	756R	700 X 480 mm



Raised rim shelf

Turned down rim

TRAYS TROLLEYS ARGON WELDED TURNED DOWN RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø100mm. With guards rails.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
435SL	600 X 400 mm	437SL	600 X 400 mm
735SL	700 X 500 mm	737SL	700 X 500 mm
835SL	800 X 580 mm	837SL	800 X 580 mm



Ref. 737SL

TRAYS TROLLEYS ARGON WELDED RAISED RIM SHELVES

Made of 18/10 stainless steel, argon welded.
Castor Ø80 mm. Delivered without accessories.

2 TRAYS		3 TRAYS	
REF.	SIZE	REF.	SIZE
351	400 X 325 mm		
356	400 X 325 mm		

OPTION : Serum rod
with clamp

REF. DESIGNATION
609I 2 stainless steel hooks

REF. DESIGNATION
609P 2 Polyamide hooks



Ref. 351

OPTIONS

Stainless steel drawers

FOR RAISED RIM SHELVES		TURN DOWN RIM SHELVES	
REF.	SIZE	REF.	SIZE
163	L240 X P370 X H85 mm	161	L240 X P370 X H85 mm

Stainless steel end-shelf

FOR TROLLEY		FOR TROLLEY	
REF.	SIZE	REF.	SIZE
4471	380 X 325 mm	7571	.450 X 325 mm

- RF** 2 castors with brake
- RAS** 4 antistatic castors Ø 100 mm
- PC** 4 ring-shaped bumpers
- 112** Holder and basin Ø 200mm
- 113** Holder with tray
- 114** Swivel plastic bag holder 10 L
- 120** Swivel plastic bag holder 30 L
- 123** Container holder



Ref. 431N with RAS et PC




Ref. 436B with RAS et PC

TRAYS TROLLEYS ARGON WELDED WITH "COMPACT" SHELVES

Made of 18/10 stainless steel, argon welded with "compact labo tacheeclean" shelves.
 Castor Ø100mm. Delivered without accessories.

2 TRAYS WITHOUT GUARD RAILS		3 TRAYS WITHOUT GUARD RAILS	
REF.	SIZE	REF.	SIZE
431	600 X 400 mm	436	600 X 400 mm
731	700 X 480 mm	736	700 X 480 mm



:
 Color option, indicate after the trolley reference the color requested :
A for Aigue marine (green)
B for bleu Oslo (blue)
N for Nectarine (orange)
P for Pamplemousse (yellow)

OPTIONS

Stainless steel end-shelf

FOR TROLLEY 600x400		FOR TROLLEY 700x48	
REF.	SIZE	REF.	SIZE
4471	L.380 X larg.325 mm	7571	L.450 X larg.325 mm

- RF** 2 castors with brake
- RAS** 4 antistatic castors Ø 100 mm
- PC** 4 ring-shaped bumpers
- 112** Holder and basin Ø 200mm
- 113** Holder with tray
- 114** Swivel plastic bag holder 10 L
- 120** Swivel plastic bag holder 30 L
- 123** Container holder

Stainless steel drawer


FOR TROLLEY WITH COMPACT TRAY	
REF.	SIZE
165	240 X 370 X 85 mm



Ref. 113



Ref. 123 (without box)



OPTIONS

FOR ALL TYPES OF TROLLEYS



■ Holder and stainless steel waste bin Réf.110



■ Gas holder Ref.111



■ Twin bags holder 10 L Ref.220



■ Holder with basin
Ø 200 mm 2 liters
Ref.112



■ Holder with tray
260x180 mm
Ref.113



■ Swivel plastic bag holder 30 L Ref. 120



■ Guard rail and clip
ref. P27
(See page 3)



■ Swivel plastic bag holder Ref.114



■ Stainless steel drawers for shaped trays
Easy to set on 600x400 or 750 500 trolleys
Two small models can be set side by side.

Small model usefull dim.255x400x115 mm
Ref.151 (for 600 x 400) Ref.153 (for 750 x 500)

Large model usefull dim.550x400x115 mm
Ref.152 (for 600 x 400) Ref.154 (for 750 x 500)



■ Drop bottles rack
For trolleys 750x500mm
Ref. 752P



■ Probe/catheter holder
Ref.508

OPTIONS FOR TROLLEYS

Distributed by
Pat Coakley Medical Ltd
sales@patcoakleymedical.ie
www.patcoakleymedical.ie
phone: 01 2018070 / 2

