



PHBOW003

General Purpose Bowl 1L

NHS Code: FSL298

✓ Pulp ✓ Recycled ✓ Disposable ✓ Single Use ✓ Multipurpose

General purpose 1 Litre bowl made of recycled and disposable pulp material with added waterproofing. High quality product designed for a multitude of functions such as carrying instruments, dressings as well as liquids and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Controlled pouring** to ease handling
- ✓ **Fits standard bowl support and dispenser rack**
- ✓ **Cover lids** for added infection control and hygiene
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full


Specification

- ✓ Capacity: 1000ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 240L x 185W x 80H (mm)
- ✓ Weight: 21g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 586mm L x 240mm W x 370mm H



Recyclable packaging



Keep away from sunlight



Manufacturer



Product conforms with applicable EC directives



Keep dry



Do not re-use



Packing

Outer: 200 pieces per case
Pallet: 24 cases per pallet



BSI Kitemark for compliance with PAS29:1999

PAS 29:1999
KM 581042

Barcodes



5 024951 012689 >

PHBOW003/APR2012

Distributed by

Pat Coakley Medical Ltd

sales@patcoakleymedical.ie

www.patcoakleymedical.ie

phone: 01 2018070 / 2

