

PHKID004

Kidney Dish

NHS Code: FTA019

✓ **Pulp** ✓ **Recycled** ✓ **Disposable** ✓ **Single Use** ✓ **Multipurpose**

Kidney Dish made of recycled and disposable pulp material with added waterproofing. High quality product designed for a multitude of functions such as carrying instruments, dressings and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full


Specification

- ✓ Capacity: 700ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 250L x 146W x 50H (mm)
- ✓ Weight: 18g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 550mm L x 250mm W x 408mm H



Recyclable packaging



Keep away from sunlight



Manufacturer



Product conforms with applicable EC directives



Keep dry



Do not re-use



Packing

Outer: 300 pieces per case
Pallet: 24 cases per pallet



PAS 29:1999
KM 581042

BSI Kitemark for compliance with PAS29:1999

Barcodes



5 024951 010999 >

PHKID004/APR2012

Distributed by
Pat Coakley Medical Ltd
sales@patcoakleymedical.ie
www.patcoakleymedical.ie
phone: 01 2018070 / 2

