



GN99

Nitrile Examination Gloves

- ✓ **Nitrile**
- ✓ **Powder Free**
- ✓ **Blue**
- ✓ **XS-XL**
- ✓ **Ambidextrous**
- ✓ **9 Newton**

Made of a strong stretchy Nitrile material, GN99 was developed to give increased performance and durability. Whilst designed for precision, GN99 is a very strong and a high quality flexible glove.

Attributes

- ✓ **Manufactured from** nitrile butadiene
- ✓ **Latex free** reducing the risk of hypersensitivity and allergic reactions
- ✓ **Low in residual chemicals** to reduce the risk of contact dermatitis associated with these chemicals
- ✓ **Flexible material** for improved feel and fit reducing fatigue during extended wear
- ✓ **Tear resistant** providing greater durability compared to latex or vinyl
- ✓ **Beaded cuff** for added strength and easier donning
- ✓ **Powder free** to reduce potential dust contamination

Regulatory compliance

- ✓ Complies with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - EN455 parts 1, 2, 3 and 4: Medical Gloves for Single use
- ✓ Manufactured in accordance with:
 - European Personal Protective Equipment Directive 89/686/EEC
 - CE Category: Cat III (0086)
 - EN420: Protective gloves: General requirements
 - EN374-2 Protective gloves: Chemicals and Micro organisms penetration
 - EN374-3 Protective gloves: Chemicals and Micro organisms permeation
- ✓ Tested to:
 - AQL 1.5
 - EN388
 - ASTM F1671-7
 - Food approval EC1935/2004
- ✓ Visually inspected for holes and minor defects
- ✓ Water tight test for pin holes in accordance with EN455-1 and ISO2859
- ✓ Typical force at break: 9 Newtons
- ✓ Shelf life of 3 years



Warning: Contains ZDBC and ZDEC chemicals, which can cause type IV skin irritation and contact dermatitis.

EN420	Level	EN374	Code
Dexterity	5	Sodium hydroxide 40%	K

Dimensions	XS	S	M	L	XL
Palm Width (mm)	≤ 80	80 ± 10	95 ± 10	110 ± 10	≥ 110
Length (mm)	Min 240	Min 240	Min 240	Min 240	Min 240
Thickness (mm) finger	0.13 ± 0.02	0.13 ± 0.02	0.13 ± 0.02	0.13 ± 0.02	0.13 ± 0.02
Thickness (mm) palm	0.10 ± 0.02	0.10 ± 0.02	0.10 ± 0.02	0.10 ± 0.02	0.10 ± 0.02
Thickness (mm) cuff	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02

EN420: glove measurements are determined by taking into account the behaviour of the glove material and intended use



Distributed by
Pat Coakley Medical Ltd
 sales@patcoakleymedical.ie
 www.patcoakleymedical.ie
 phone: 01 2018070 / 2






Packaging

Outer



Inner



Dimensions

Inner: 235mm L x 125mm W x 70mm H

Outer: 360mm L x 255mm W x 245mm H

Size Colour Guide



Barcodes	Size	Inner Barcode	Outer Barcode
 5 024951 014690 >	XS	5024951014683	5024951014737
	S	5024951014690	5024951014744
	M	5024951014706	5024951014751
	L	5024951014713	5024951014768
	XL	5024951014720	5024951014775



Packing

Inner: 100 gloves per dispenser

Outer: 10 dispensers per case

Pallet: 66 cases per pallet

Symbols



Lot number



Compliant with the Medical Devices Directive



Compliant with the Personal Protective Equipment Directive



Temperature limitation



Date of manufacture



Non-sterile



EN420: Protective gloves: General requirements



Recyclable packaging



Use by



Acceptable quality level



EN374 - 2: Micro organism hazards



Keep dry



Manufacturer



Latex free



EN374 - 3: Chemical hazards



Keep away from direct sunlight and fluorescent lighting



Do not re-use



Powder free



Food approved

GN99/NOV2013

Distributed by
 Pat Coakley Medical Ltd
 sales@patcoakleymedical.ie
 www.patcoakleymedical.ie
 phone: 01 2018070 / 2

