

HandSafe® GS690

Sterile Nitrile Examination Gloves

Powder Free Nitrile Examination Gloves, Sterile

Gloves Specification:

- Nitrile
- Textured
- Blue
- Beaded Cuff For added strength and easier donning
- Sterile
- Strong, soft and flexible
- Latex Free
- Ambidextrous

Packing:

Pouch: 1 pair per pouch

Dispenser: 50 pairs per dispenser box

Outer Case: 4 dispensers per case

Pallet: 54 cases per pallet

Performance Specification:

Freedom From holes - BS EN 455-1	
Water Tight Test	General Inspection Level 1, AQL 1.5

Dimensions – BS EN 455-2	Small	Medium	Large
Palm Width (mm)	80 ± 10	95 ± 10	110 ± 10
Length (mm)	240 min	240 min	240 min

Force at Break – BS EN 455-2	
Throughout shelf life	Minimum 6 Newton
After challenge test	Minimum 6 Newton

Biological Properties – BS EN 455-3	
Powder Level	Less than 2 mg per glove
Protein Content (Latex only)	≤50µg/g using the modified Lowry method

Viral Penetration – ASTM F1671-7	
Phi-X174 Bacteriophage penetration	No visual penetration

Shelf Life Properties – BS EN 455-4	
Currently shelf life is determined through accelerated ageing and is set at 3 years (whilst real time studies are in progress)	

Sterilisation Method – BS EN ISO 11137	
Sterilisation	Gamma Irradiation

Quality Standards:

Manufacturing QMS	EN ISO 9001 / EN ISO 13485
Glove	The BS EN 455 series of standards (parts 1,2,3 and 4)
Sampling & Inspection	ISO2859-1
MDD Requirement	Medical Devices Directive 93/42/EEC (and subsequent amendments)
CE Classification	Class I Sterile Medical Device
Glove Care	Keep away from direct sunlight or Fluorescent light. Store at temperature no less than 5°C, no more than 25°C
Sterility Care	Sterile unless package is opened or damaged




HandSafe® GS690

Sterile Nitrile Examination Gloves

Packaging Artwork:

Inner



Outer



Pouch:



Dimensions:

Dispenser: L145mm x W135mm x H235mm

Outer: L305mm x W278mm x H247mm

Pouch: W116mm x H116mm


Barcodes:

Size	Small	Medium	Large
Inner	5052755002209	5052755002216	5052755002223
Outer	5052755002230	5052755002247	5052755002254

Packaging Symbols:



Lot number



Compliant with the Medical Devices Directive



Keep dry



Date of manufacture



Single use - do not re-use




Keep away from sunlight



Use by



Latex Free



Temperature limitation



Manufacturer



Powder Free



Recyclable packaging



Acceptable quality level



Sterilization by Gamma Irradiation

GS690/OCT2013

Distributed by

Pat Coakley Medical Ltd

sales@patcoakleymedical.ie

www.patcoakleymedical.ie

phone: 01 2018070 / 2

