

# TRACHEOSTOMY CARE

## ADVAZORB FIXATION DRESSING



LATEX  
FREE




Advazorb Fixation is a post decannulation tracheostomy dressing fitted with a central location button that prevents unwanted escape of air through the wound site with gentle pressure applied during occasions such as talking, coughing and sneezing.

Advazorb is a hydrophilic foam dressing with a unique perforated soft silicone wound contact layer and surrounding adhesive. It is an atraumatic foam dressing which does not adhere to moist wounds.

The silicone helps maintain wound healing environment and makes it ideal as a non-adherent wound contact dressing. Soft silicone dressings do not damage newly formed granulating tissue and do not lose their adhesion after initial application.

The silicone border prevents the need for additional dressing or fixation. Silicone adhesive combines security with pain-free dressing changes. Waterproof backing enables the patient to shower without the need to replace or cover the dressing.

- Supplied sterile and individually wrapped, ready for use
- Location button enables accurate identification of tracheostomy site
- Maintains a moist environment and manages exudates
- Allows atraumatic dressing change
- Provides a barrier for bacteria
- Latex-Free
- UK Manufactured
- CE Marked



IMP\$CODE	DESCRIPTION	SIZE	DRUG\$TARIFF	NHS\$CODE	PACK\$QTY
CR/4341	Advazorb/Fixation	12cm/x/10cm	9	FDH605	10