

# TRACHEOSTOMY CARE

## ADVAZORB PLUS (T) DRESSING



LATEX  
FREE




Advazorb Plus (T) offers a new standard in hydrophilic foam tracheostomy dressings with incredible absorbency ability.

Manufactured in lightweight, soft and conforming low-adherent foam with high fluid absorbance and retention characteristics.

Newly improved with a film backing added to ensure improved strength and absorbance for medium/high exuding wounds.

Supplied sterile and individually wrapped, ready for use. Available in one convenient size to suit most intubation procedures, and to offer comfort to the patient.

- Cushioning effect on wound
- Low adherence -Reduced trauma on removal
- Highly absorbent 5mm foam
- Rapid uptake of exudates
- Can be cut to size
- Latex-Free
- UK Manufactured
- CE Marked



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3900	Advazorb Plus Dressing	80mm x 80mm	-	FDH419	20



**Pat Coakley Medical**  
E3 & E4 Network Business Park  
Kilcoole, Co Wicklow

Tel: +353 1 201 8070/2/3  
Fax: +353 1 201 8074  
[www.patcoakleymedical.ie](http://www.patcoakleymedical.ie)