

TUBE FIXATION

TUBE TIDY CLIPS




The simplest product on the market for keeping tubes and lines arranged and sorted during patient monitoring.

The product can clip to blankets, sheets, pillows and opens to allow tubing to run through.

Simply lay tubes required over hook/loop and seal to hold tubes in place.

Organises cords and drapes, useful in all specialist units and wards.

- Can retain drapes and linens without the use of pins or fixers
- User is able to write on the product for identity markings
- Simply open the hook/loop material, insert tubing and close for use
- Latex-Free
- UK Manufactured
- CE Marked



IMP CODE	DESCRIPTION	COLOUR	DRUG TARIFF	NHS CODE	PACK QTY
TFP/800	Tube Tidy Clips	White	-	-	10
TFP/850.B	Tube Tidy Clips	Blue	-	-	20
TFP/850.R	Tube Tidy Clips	Red	-	-	20
TFP/850.Y	Tube Tidy Clips	Yellow	-	-	20



Insight Medical Products Ltd
Unit 3, Priory Industrial Estate,
Tetbury, Gloucestershire. GL8 8HZ

Tel: +44 (0)1666 500055
info@insightmedical.net
www.insightmedical.net

