



COLLECTION



# CALZURO® Classic

---



Professional clogs, made of a thermo-tecno polymer compound latex-free, completely antistatic, with side holes for foot perspiration, available with or without upper holes, autoclavable to 134°C/273°F, washable to 90°C/194°F.

Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20347:2012, requirements:**

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRA** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1314685**

## SIZES

34/35 – 35/36 – 36/37 – 37/38 – 38/39 – 39/40 – 40/41 – 42/43 – 44/45 – 46/47

## COLORS



WHITE



GREEN



LIGHT BLUE



ORANGE



BLACK



PINK <sup>(2)</sup>



DARK BLUE



MET. BLUE



RED



PURPLE



MAROON



YELLOW



MET. GREY



TEAL <sup>(1)</sup>

(1) available only with upper holes

(2) available only up to 42/43

AVAILABLE ALSO WITH HEEL-STRAP



CE CERTIFIED



ECO-FRIENDLY



100% RECYCLABLE



LATEX-FREE



ANTISLIP



AUTOCCLAVABLE



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 90°  
WASHING MACHINE



UNISEX

# ACCESSORIES **Classic**

---

## INSOLE



Arch support insole for Calzuro, in special foam that allows feet to breathe, autoclavable to 121°C/250°F up to 50 times.

## SIZES

34 – 35 – 36 – 37 – 38 – 39 – 40 – 42 – 44 – 46

## CHARMS




## CALZURELLI



## SWAROVSKI



HEEL STRAP  
ON DEMAND



COLORS



**NB:** FANCY COLORS SUPPLIED ONLY WITH WHITE HEEL STRAP  
Autoclavable to 134°C/273°F

CALZURO® **Fancy**

---



Professional clogs, made of a thermo-tecno polymer compound latex-free, completely antistatic, with side holes for foot perspiration, available only with upper holes, autoclavable to 134°C/273°F, washable to 90°C/194°F.

Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20347:2012, requirements:**

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRA** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1314685**

## SIZES

34/35 – 35/36 – 36/37 – 37/38 – 38/39 – 39/40 – 40/41 – 42/43 – 44/45 – 46/47

## COLORS



FANCY  
WHITE-BLUE<sup>(1)</sup>



FANCY  
WHITE-PINK<sup>(1x2)</sup>

(1) available only with upper holes

(2) available only up to 42/43

SUPPLIED ONLY WITH WHITE HEEL-STRAP



CE CERTIFIED



ECO-FRIENDLY



100% RECYCLABLE



LATEX-FREE



ANTISLIP



AUTOC LAVABLE



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 90°  
WASHING MACHINE



UNISEX

CALZURO® **Light**

---





Professional clogs, made of ultra-light soft plastic latex-free material, casein-resistant, with antistatic plug and side holes for foot perspiration. Available only without upper holes and with back heel-strap. Machine washable up to 50°C/122°F. Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20347:2012, requirements:**

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1314687**

## SIZES

35 – 36 – 37 – 38 – 39 – 40 – 41 – 42 – 43 – 44 – 45 – 46

## COLORS



WHITE



GREEN



NAVY BLUE



BLACK



PINK

(only up to size 41) (only up to size 44)



MAUVE

NEW



LIGHT BLUE

NEW



OLIVE GREEN



CE CERTIFIED



ECO-FRIENDLY



LATEX-FREE



ANTISLIP



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 50°  
WASHING MACHINE



UNISEX

# ACCESSORIES **Light**



## ANTISTATIC INSOLE

Comfortable ANTISTATIC insole, Drytech® fabric keeps the foot dry and latex foam bottom provides high level of comfort. Antibacterial treatment can prevent bacteria causing unpleasant odors.

### SIZES

35 – 36 – 37 – 38 – 39 – 40  
41 – 42 – 43 – 44 – 45 – 46



## CUSTOM FIT INSOLE

Comfortable anti-perspiring insole, which **can be trimmed for custom fit**.

Drytech® fabric keeps the foot dry and latex foam bottom provides high level of comfort.

Antibacterial treatment can prevent bacteria causing unpleasant odors.

Machine washable up to 30°C/86°F




## PILE INSOLE

Comfortable and warm PILE insole, which **can be trimmed for custom fit**.

Drytech® fabric keeps the foot dry and latex foam bottom provides high level of comfort. Antibacterial treatment can prevent bacteria causing unpleasant odors.

Machine washable up to 30°C/86°F

## HEEL-STRAP REPLACEMENT



## COLORS



WHITE



GREEN



NAVY BLUE



ORANGE



BLACK



PINK



MAUVE

# CALZURO® **Pro Safety**

---

(WITH SAFETY TOE-CAP)



Soft and practical professional footwear with closed back-heel, made of a brand-new special ultra-light copolymer, latex-free, casein-resistant.

Easy wear slip on and slip off and safe fit for maximum comfort and further your safety policy at work. Supplied without side vents, with antistatic plug and with synthetic polymer protection metal-free toe-cap (*impact resistance of the toe-cap: 200 Joule on the toe*).

The safety toe-cap provides adequate protection against mechanical risks, such as impact and crushing/breakthrough, but it shows greater lightness than the corresponding leather shoes with steel toe-cap.

Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20345:2011, requirements:**

**S4** (additional requirements):

- **Heel area closed**
- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **C leated outsole**

**SRC** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1770379**

SIZES

36 – 37 – 38 – 39 – 40 – 41 – 42 – 43 – 44 – 45 – 46

COLORS



CE CERTIFIED



ECO-FRIENDLY



ANTISLIP



LATEX-FREE



SAFETY TOE-CAP  
METAL-FREE



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 50°  
WASHING MACHINE



UNISEX



EASY-FIT

CALZURO® **Pro SP**

---

(WITHOUT SAFETY TOE-CAP)



Soft and practical professional footwear with closed back-heel, made of a brand-new special ultra-light copolymer, latex-free, casein-resistant.

Easy wear slip on and slip off and safe fit for maximum comfort and further your safety policy at work.

Supplied with side vents for foot perspiration and with antistatic plug.

Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20347:2012, requirements:**

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1639834**

## SIZES

36 – 37 – 38 – 39 – 40 – 41 – 42 – 43 – 44 – 45 – 46

## COLORS



WHITE



NAVY BLUE



RASPBERRY  
(ONLY UP TO SIZE 42)



CE CERTIFIED



ECO-FRIENDLY



ANTISLIP



LATEX-FREE



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 50°  
WASHING MACHINE



UNISEX



EASY-FIT

# CALZURO® Aqua

Flip flops made of EVARITE®, a new, soft, lightweight, latex-free plastic material which does not favor the proliferation of fungi and bacteria.

Anatomically shaped footbed, well-dressing anatomic shape, thong comfortable to wear, slightly raised heel and toe bed.

Machine washable up to 50°C/122°F. Slip-resistant sole CE certified. 100% Made in Italy.

**Certified according to EN ISO 20347:2012**, requirements:

- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

**MEDICAL DEVICE CLASS I NON STERILE: DM 1587673**

## SIZES

35 – 36 – 37 – 38 – 39 – 40  
41 – 42 – 43 – 44 – 45 – 46

## COLORS



BLACK



TURQUOISE



NAVY BLUE



BRIGHT GREEN



AMETHYST



FUXIA

(ONLY UP TO SIZE 41) (ONLY UP TO SIZE 41)



CE CERTIFIED



ECO-FRIENDLY



LATEX-FREE



ANTISLIP



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 50°  
WASHING MACHINE



UNISEX



# CALZURO® Buddy

Latex-free, antistatic slippers, washable up to 90°C/194°F, autoclavable to 121°C/250°F, with upper and side holes for foot perspiration, insulating sole from cold and hot, anatomical inner arch-support. Supplied with adjustable heel-strap. 100% Made in Italy.

## SIZES

35 - 36 - 37 - 38 - 39 - 40/41 - 42/43 - 44/45

## COLORS



WHITE +  
METALLIC  
PURPLE\*

\*(only up to  
size 42/43)



WHITE + BLACK



WHITE + METALLIC  
BLUE



100% RECYCLABLE



LATEX-FREE



AUTOC LAVABLE



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



WASHABLE IN 90°  
WASHING MACHINE



UNISEX

# CALZURO® Bab

Antistatic white boots in special polymer compound (*nitrile rubber*), non toxic, resistant to hydrocarbons (*fuel oil*). A good washing resistance at medium-high temperatures 85°C/185°F helps maintain within optimal limits to the spread of bacteria. Resistant to animal fats, vegetable oils and casein, slip-resistant sole. Available: model at knee.

White polyester knitting inside.

Washable up to 85°C/185°F. 100% Made in Italy.

**BAB** (to the knee)

**Certified according to EN ISO 20347:2012**, requirements:

- **O4** (closed seat region, antistatic properties, energy absorption of seat region)
- **FO** (bottom sole resistance to fuel oil)
- **SRC** (Slip resistance)

## SIZES

36 – 37 – 38 – 39 – 40 – 41 – 42  
43 – 44 – 45 – 46 – 47 – 48

## COLORS



WHITE



CE CERTIFIED



ANTISLIP



FUEL OIL - RESISTANT



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



WASHABLE IN 85°  
WASHING MACHINE



MADE IN ITALY



UNISEX



White boots in special polymer compound (*nitrile rubber*), non toxic, model at calf, white polyester kitting inside.  
Washable up to 85°C/185°F. 100% Made in Italy.



### **BAB Half Wellington (to the calf)**

**Certified according to EN ISO 20347:2012**, requirements:

- **OB** (basic requirements)
- **SRA** (slip resistance)

### **SIZES**

36 – 37 – 38 – 39 – 40 – 41  
42 – 43 – 44 – 45 – 46 – 47

### **COLORS**



WHITE



CE CERTIFIED



MADE IN ITALY



UNISEX



ANTISLIP



FUEL OIL - RESISTANT



CRUELTY-FREE



WASHABLE IN 85°  
WASHING MACHINE

# CALZURO® Effetto

Slippers for medical and paramedical Personnel, slip-resistant PU (Polyurethane) sole, leather upper made of pigmented perforated leather for a healthy and natural foot perspiration, soft padding on the instep. Can be supplied with back heel-strap.

**Certified according to EN ISO 20347:2012**, requirements:

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

## MODELS



WITHOUT HEEL-STRAP



WITH HEEL-STRAP



## SIZES

35 – 36 – 37 – 38 – 39 – 40 – 41  
42 – 43 – 44 – 45 – 46

## COLORS



WHITE  
(white sole)



BLUE  
(black sole)



ANTISTATIC



UNISEX



MADE IN ITALY



CE CERTIFIED



ANTISLIP

# CALZURO® Giulia

Women's leather slippers with removable double-face insole\* (massaging/leather) and upper heel-strap to adjust the instep, TR (thermoplastic material) bottom sole.

## MODELS



WITH UPPER HOLES  
AND WITHOUT HEEL-STRAP



WITH UPPER HOLES  
AND WITHOUT HEEL-STRAP



WITHOUT UPPER HOLES  
AND WITH HEEL-STRAP



\*ALWAYS SUPPLIED WITH  
FOOTWEAR ONLY

## SIZES

35 - 36 - 37 - 38 - 39 - 40 - 41

## COLORS



WHITE  
(white sole)



BLUE  
(Available only with upper holes  
and without heel-strap)



# CALZURO® **Wood**

Calzuro® Wood the “flexible”, with perforated upper in pure Italian leather and anti-bacterial fabric inside, flexible sole in beech wood and special rubber.

Washable up to 40°C/104°F.

## SIZES

35 - 36 - 37 - 38 - 39 - 40  
41 - 42 - 43 - 44 - 45 - 46

## COLORS



WHITE



NAVY BLUE



WASHABLE IN 40°  
WASHING MACHINE



UNISEX

# CALZURO® **Comfort**

---

Comfort Slippers entirely made of leather and sewed by hand, with back heel-strap, anatomical insole in natural perspiring leather, TR (thermoplastic material) bottom sole.

## SIZES

---

35 - 36 - 37 - 38 - 39 - 40 - 41

## COLORS

---



WHITE



BLUE





**100% MADE IN ITALY**

[www.calzuro.it](http://www.calzuro.it)

**Calzuro®** is a registered trademark owned by BIHOS Srl

[www.bihos.com](http://www.bihos.com)



*Sole Distributor in Ireland*

**Pat Coakley Medical Ltd**

Units E3 E4  
Network Enterprise Park  
Kilcoole A63 H244  
Co Wicklow  
IRELAND

+353 1 2018070  
[sales@patcoakleymedical.ie](mailto:sales@patcoakleymedical.ie)