

GN90

Blue nitrile powder free disposable glove

GN90



Sectors:

- Automotive
- Food Manufacturing and Preparation
- Janitorial
- Medical, Vet and Dental

Features

Manufactured from: Specially formulated nitrile butadiene polymer

Latex free: reducing the risk of hypersensitivity and allergic reactions

Low in residual chemicals: To reduce the risk of contact dermatitis associated with these chemicals

Flexible material: For improved feel and fit reducing fatigue during extended wear

Strength: Offers better puncture, tear and abrasion resistance than natural rubber or vinyl

Food Contact: Tested in accordance with European standard EN1186 making it suitable for food contact

Powder free: Reduces the risk of dust contamination

Textured finish: For improved grip

Pat Coakley Medical 



MDDCLASS I
PPE CAT III
2777



K



VIRUS



Code	GN90
Colour	Blue
Sizes	XS-XL



Inners dimensions (mm)
235 x 125 x 82mm

Outers dimensions (mm)
418 x 258 x 251mm

Size	Inner Barcode	Outer Barcode
XXS/6	5024951012719	5024951012764
S/7	5024951012726	5024951012771
M/8	5024951012733	5024951012788
L/9	5024951012740	5024951012795
XL/10	5024951012757	5024951012801

Standards

Contains accelerators which may cause allergic reactions

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with the Medical Devices Directive 93/42/EEC Class 1

Tested to European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Product Food Approved to EC1935/2004

EU Type Examination carried out by SATRA

Tested to European Standards EN420:2003+A1:2009 & EN ISO 374-1:2016 & EN ISO 374-5:2016

Pat Coakley Medical Ltd.

Units E3 E4 Network Enterprise Park
Kilcoole A63 H244 Co Wicklow
IRELAND

+353 12018070
sales@patcoakleymedical.ie

Packing

200 gloves per dispenser
10 dispensers per case
60 cases per pallet

For further technical information, please contact Customer Services on 033 3320 8550 or the Technical team on 033 3320 8555

