



Vitrex® Blood gas capillary Tubes

Vitrex disposable capillary tubes for blood gas analyses are made from heavy walled glass, suitable for 2.3 mm Ø adaptors.

Uniform heparinisation, careful fire polishing, heavy walled tubes and robust end-caps secure proper function and reduce operating risk. The strict demands on function and quality for the VITREX products are your guarantee.



Pat Coakley Medical Ltd
 Units E3 E4, Network Enterprise Park, Kilcoole, A63 H244 Co Wicklow, IRELAND
 Phone +353 12018070 sales@patcoakleymedical.ie

Packaging

Choose your Vitrex blood-gas-capillary tubes according to diagram shown below. Just indicate article number and type.

The types have been manufactured under strict control and are in accordance with current standards available in 3 lengths: 75 mm, 100 mm and 125 mm. Packed in vials of 250 pcs. Special types upon request.

Accessories

Stirrers packed in vials of 250 pcs.
 End-caps packed in boxes of 500 pcs.



Technical Specs.

Vials of	Na Hep. Article No.	Li Hep. Article No.	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 1 mm	Vol. (µl)
4 x 250	182513	182613	1.6	2.3	75	130 (min)
4 x 250	182413	183313	1.6	2.3	100	175 (min)
4 x 250	182313	182013	1.6	2.3	125	220 (min)

Shelf life 4 years

Accessories	Article No.	Description
10 x 500	140813	End caps
10 x 250	140913	Stirrers
1	142400	Magnets

