

Trans-anal Irrigation: Solutions with you in mind



Rectal Irrigation or Trans-anal Irrigation (TAI)

Effective bowel emptying is the norm for most people. For some, however, functional or neurogenic bowel disorders may cause faecal incontinence and chronic constipation. These symptoms are often distressing and can have a negative effect on a person's lifestyle.

Regular bowel emptying through use of Rectal Irrigation, also known as Trans-anal Irrigation (TAI), can help people avoid complications and maintain their quality of life.






Renew Aquaflush[®] Irrigation system

Solutions with you in mind

All **Renew Aquaflush** Irrigation systems are modular, resulting in a solution that can be modified according to individual need. Our Clinical Advisory team are skilled to work in partnership with healthcare professionals and patients to determine the most suitable system.

Lite

Aquaflush Lite is our simple full irrigation system that uses modular components to deliver an intuitive and easy-to-use solution. Featuring a large 1.2 litre bag, it is suitable for patients who can sit on a toilet or commode.



Features

- 1 Super-smooth silicone cone with rounded tip, side eyelets & ergonomic grip
- 2 Washbag for discreet storage at home & when travelling
- 3 1.2-litre water bag with screw-cap & thermometer ensures safe water temperature
- 4 Hand pump to easily control speed of irrigation
- 5 In-line tap to control water flow
- 6 'Click/close' connector with one-way valve prevents dirty water entering system



Actif

Aquaflush Actif is a unique dual-purpose product for either paediatric rectal irrigation or adult stoma irrigation. The product features an ergonomically designed cone which is the perfect size to suit both applications.



Compact & New Compact +

Aquaflush Compact & Compact + have been designed with you in mind by providing a system that is quick, easy to learn and safe to use. Compact is designed for those who only need a small amount of water when irrigating.

Accessories



Water-bag stand



Extension tube

Accessories (continued)



Lubricating jelly



Cone Refill Pack



Short Cone Refill Pack

RENEW Aquaflush order codes

Aquaflush Product Description	Product Code	GMS Code
Compact Starter Kit	AFCS	93967
Compact Monthly Kit	AFCM	93966
Compact + Starter Kit	AFCPS	95003
Compact + Monthly Kit	AFCPM	95002
Lite Starter Kit	AFLS	93972
Lite Monthly Kit	AFLM	93968
Actif Starter Kit	AFAS	94229
Actif Monthly Kit	AFAM	93970

Aquaflush Accessories Description		
Lite Cones x 15	AFLA	94357
Short Cones x 15	AFSC	94945
Extension Tubes x 2	AFEXTN	93969

Pat Coakley Medical 

For more information on **AQUAFLUSH**,
visit the RENEW Medical website at
www.renew-medical.uk or call
Pat Coakley Medical on **01 2018070**
or **086 7741247**

 **Renew**

Renew Medical Limited
15 Southampton Place
London, WC1A 2AJ
Freephone: 0800 542 0814
Email: support@renew-medical.uk
web: www.renew-medical.uk