


PURIFOAM LARYNGECTOMY FILTERS

Purifoam Filters use a high quality medical grade foam that protects the stoma from dirt and dust entering the airway, while warming and humidifying the incoming air.

Purifoam is adhered to the skin using a skin friendly adhesive strip that has rounded edges for added comfort.

Manufactured from beige coloured foam to give a natural aesthetic look. Purifoam has excellent moisture retaining qualities, ensuring good levels of absorbency to give exceptional humidification to the airway for breathing reassurance.



- Change daily or more frequently if required
- Available on FP10 Prescription
- Single patient use only
- UK Manufactured
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3853	Purifoam Filters	6.5cm x 6.5cm	319-0386	FTH029	30