

SOFSHIELD LARYNGECTOMY & TRACHEOSTOMY PROTECTORS



Sofshield Laryngectomy & Tracheostomy Protectors can be washed a minimum of 10 times, while still retaining filtration efficiency.

This offers a significant cost saving and longer product use over alternative protectors that use foam inserts, while also eliminating the need for costly additional cosmetic covers.

All Sofshield Protectors have an antibacterial coating to help prevent infections.

Sofshield is available in five colours and three sizes to suit all Laryngectomy and Tracheostomy patients.

Sofshield uses a silicone toggle system that eliminates the need to tie cords, therefore allowing patients with restrictive movement an easier fitting process.



- Antibacterial coating to prevent infection
- Consistent filtration efficiency after being washed
- Available on FP10 Prescription
- Single patient use only
- UK Manufactured
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
AS/3669	Sofshield -White	Large	291-9231	FTH006	10
AS/3764	Sofshield -Blue	Large	305-5373	FTH008	10
AS/3765	Sofshield -Beige	Large	305-5381	FTH007	10
AS/3698	Sofshield -White	Small	291-9249	FTH009	10
AS/3777	Sofshield -Blue	Small	305-5399	FTH100	10
AS/3778	Sofshield -Beige	Small	305-5407	FTH010	10
AS/4051	Sofshield -Black	Medium	377-1722	FSM2621	10
AS/4056	Sofshield -Light Blue	Medium	409-1526	FTH129	10