


LARYNGECTOMY CARE

VOICE PROSTHESIS BRUSHES

LATEX FREE  

Voice Prosthesis Brushes are manufactured with soft nylon bristles that prevent damaging the silicone inner surface of the prosthesis.

Brushes are tapered to provide effective cleaning of mucus and secretion build-up in the inner lumen of the voice prosthesis.



- Silicone safety sleeve to prevent over-insertion
- Flexible plastic handle for ease of use
- Listed on FP10 Drug Tariff
- Latex-Free
- UK Manufactured
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
PB/50	Voice Prosthesis Brush	Standard	22974395	7	10