



# PHBED002

## Bedpan Liner

NHS Code: FSE063

✓Pulp ✓Recycled ✓Disposable ✓Single Use

Bedpan liner made of recycled and disposable pulp material with added waterproofing. High quality product designed for enhanced patient acceptance and comfort. Designed to be used with a bedpan holder and certain commode chairs, suitable for disposal in any macerator.

### Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Fits standard supports, holders and dispensing racks**
- ✓ **Controlled pouring** to ease handling
- ✓ **Cover lids** for added infection control and hygiene
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full

### Specification

- ✓ Capacity: 2000ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 355L x 296W x 95H (mm)
- ✓ Weight: 52g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
  - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
  - British standard PAS29:1999
- ✓ Disposability: By maceration

### Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

### Packaging



Dimensions: 368mm L x 305mm W x 623mm H

- Recyclable packaging
- Keep away from sunlight
- Manufacturer
- Product conforms with applicable EC directives
- Keep dry
- Do not re-use



### Packing

Outer: 100 pieces per case  
Pallet: 18 cases per pallet



PAS 29:1999  
KM 581042

BSI Kitemark for compliance with PAS29:1999

### Barcodes



PHBED002/APR2012

Distributed by  
**Pat Coakley Medical Ltd**  
sales@patcoakleymedical.ie  
www.patcoakleymedical.ie  
phone: 01 2018070 / 2

