



PHURI046

Bulbous Male Urinal

NHS Code: FWV042

✓Pulp ✓Recycled ✓Disposable ✓Single Use

Bulbous male urinal made of recycled and disposable pulp material with added waterproofing. High quality product for enhanced patient acceptance and comfort. Designed to be used with a bedside urinal holder and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Controlled pouring** to ease handling
- ✓ **Fits** standard dispensers and holders
- ✓ **Smooth finished opening** to enable easy pouring and for ergonomic purposes
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Easy to carry** when full


Specification

- ✓ Capacity: 875ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 260L x 105W x 125H (mm)
- ✓ Weight: 38g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 842mm L x 252mm W x 900mm H



Recyclable packaging



Keep away from sunlight



Manufacturer



Product conforms with applicable EC directives



Keep dry



Do not re-use



Packing

Outer: 100 pieces per case
Pallet: 5 cases per pallet



BSI Kitemark for compliance with PAS29:1999

PAS 29:1999
KM 581042

Barcodes



5 024951 012689 >

PHURI046/APR2012

Distributed by

Pat Coakley Medical Ltd

sales@patcoakleymedical.ie

www.patcoakleymedical.ie

phone: 01 2018070 / 2

