



# PHPOT063

## Waxed Paper Medicine Pot

NHS Code: FTP087

✓ Paper ✓ Recycled ✓ Disposable ✓ Single Use ✓ Multipurpose

Medicine pot made of disposable paper with added waxing for light waterproofing. High quality product designed for a multitude of functions but primarily used for dispensing pills and suitable for disposal in nay macerator.

### Benefits

- ✓ Disposable product for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ Pulp product made of recycled materials and disposable through macerators
- ✓ Fluid Resistant to allow use with liquids
- ✓ Steady when standing due to flat bottom surface
- ✓ Reinforced rolled edges for increased rigidity
- ✓ Easy to store (product is stackable) and separate
- ✓ Easy to carry when full

### Specification


- ✓ Capacity: 28ml
- ✓ Material: Waxed paper
- ✓ Dimensions: 40L x 40W x 29H (mm)
- ✓ Weight: 0.1g ± 4g
- ✓ Colour: White
- ✓ Disposability: By maceration

### Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)



### Packaging



Dimensions: 89mm L x 45mm W x 343mm H



Dimensions: 359mm L x 229mm W x 345mm H

CE Product conforms with applicable EC directives

Recyclable packaging

Keep away from sunlight

Do not re-use

Manufacturer

Keep dry

### Packing



Inner: 250 pots per inner  
Outer: 5000 pots per case  
Pallet: 30 cases per pallet

### Barcodes



Inner	5024951012481
Outer	5024951012498

PHPHOT063/APR2012



0044 (0) 208 8335 3636

0044 (0) 208 8335 0300

@ sales@hpchealthline.co.uk

HPC Healthline Limited, Colwood House, 158 Garth Road, Morden, Surrey, SM4 4LZ, England

www.hpchealthline.co.uk

www.facebook.com/hpchealthline

#hpchealthline