



PHCUP008

Multi Purpose Cup

NHS Code: FSX016

✓Pulp ✓Recycled ✓Disposable ✓Single Use ✓Multipurpose

Multi purpose cup made of recycled and disposable pulp material with added waterproofing. High quality product designed for a multitude of functions, but primarily used for urine samples and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Controlled pouring** to ease handling
- ✓ **Cover lids** for added infection control and hygiene
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full


Specification

- ✓ Capacity: 300ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 176L x 90W x 64H (mm)
- ✓ Weight: 10g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 392mm L x 176mm W x 431mm H



Recyclable packaging



Keep away from sunlight



Manufacturer



Product conforms with applicable EC directives



Keep dry



Do not re-use



Packing

Outer: 200 pieces per case
Pallet: 36 cases per pallet



BSI Kitemark for compliance with PAS29:1999

PAS 29:1999
KM 581042

Barcodes



5 024951 012856 >

PHCUP008/APR2012



0044 (0) 208 8335 3636



0044 (0) 208 8335 0300



sales@hpchealthline.co.uk



HPC Healthline Limited, Colwood House, 158 Garth Road, Morden, Surrey, SM4 4LZ, England



www.hpchealthline.co.uk



www.facebook.com/hpchealthline



#hpchealthline

HPC Healthline