




## ADVAZORB (T) LITE TRACHEOSTOMY DRESSING

Advazorb (T) Lite is the definitive lightweight hydrophilic slim foam tracheostomy dressing.

Advazorb (T) offers remarkable standards of conformability, softness, lightness compared to gauze dressings.

Perfect for the treatment of light/medium exuding wounds as moisture is absorbed and retained in the hydrophilic cellular structure.

Supplied sterile and individually wrapped, ready for use. Available in one convenient size to suit most intubation procedures, and to offer comfort to the patient.



- Cushioning effect on wound
- Low adherence -Reduced trauma on removal
- Highly absorbent 2.5mm foam
- Excellent patient comfort
- Rapid uptake of exudates
- Available on FP10 Prescription
- Single patient use only
- UK Manufactured
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/4212	Advazorb Lite Dressing	80 x 80mm	374-0230	ELY392	20