

## ADVAZORB PLUS (T) TRACHEOSTOMY DRESSING




Advazorb Plus (T) offers a good standard in hydrophilic foam tracheostomy dressings with high absorbency ability.

Manufactured in lightweight, soft and conforming low-adherent foam with high fluid absorbance and retention characteristics.

Manufactured with a film backing to ensure improved strength and absorbance for medium/high exuding wounds.

Supplied sterile and individually wrapped, ready for use.

Available in one convenient size to suit most intubation procedures, and to offer comfort to the patient.



- Cushioning effect on wound
- Low adherence -Reduced trauma on removal
- Highly absorbent 5mm foam
- Excellent patient comfort
- Rapid uptake of exudates
- Can be cut easily to size
- Available on FP10 Prescription
- Single patient use only
- UK Manufactured
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3900	Advazorb Plus Dressing	80 x 80mm	395-1126	FDH419	20