

LARYNGECTOMY CARE

PURIFOAM LARYNGECTOMY FILTERS



LATEX
FREE




Purifoam Filters use a high quality medical grade foam that protects the stoma from dirt and dust entering the airway, while warming and humidifying the incoming air.

Purifoams are adhered to the skin using a skin friendly adhesive strip that has rounded edges for added comfort.

Manufactured from beige coloured foam to give a natural aesthetic look. Purifoam has excellent moisture retaining qualities, ensuring high levels of absorbency to give exceptional humidification to the airway for breathing reassurance.

- Change daily or more frequent if required
- Listed on FP10 Drug Tariff
- Latex-Free
- UK Manufactured
- CE Marked



IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
CR/3853	Purifoam Filters	6.5cm x 6.5cm	319-0386	FTH029	30



Pat Coakley Medical
E3 & E4 Network Business Park
Kilcoole, Co Wicklow

Tel: +353 1 201 8070/2/3
Fax: +353 1 201 8074
www.patcoakleymedical.ie