

LARYNGECTOMY PROTECTORS

CASCADE SHOWER PROTECTOR



LATEX
FREE



Cascade is a soft waterproof PVC protector that provides a barrier to prevent shower water from entering a stoma thereby improving the patient's confidence to shower routinely.

The neckband is soft and has elastic to form a comfortable yet firm seal around the neck to prevent water seeping below the protector. The backing is a knitted polyamide material that will absorb moisture from surrounding wet skin and provide a comfortable warm feel, unlike other products that may initially feel cold to the skin.

From a safety perspective, the polyamide backing also prevents adhesion to the skin or stoma button when bib is wet, as adhesion from unlined protectors can reduce the ability for a patient to breath properly.

For added comfort and breathability the Cascade Protector can be worn over a Sofshield Protector or alternative filter, however, this is not essential due to the backing material. The protector is an aesthetically pleasing white in colour and supplied in pairs to allow for drying between showers.



- Shower Proof
- Lay flat contour
- Washable
- Listed on FP10 Drug Tariff
- Latex-Free
- UK Manufactured
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
AS/3835	Cascade Shower Protector	Standard	314-2601	FTH028	10



Pat Coakley Medical
E3 & E4 Network Business Park
Kilcoole, Co Wicklow

Tel: +353 1 201 8070/2/3
Fax: +353 1 201 8074
www.patcoakleymedical.ie