

Plastic Holloware

Cape Warwick

Cape Warwick Ltd. are the original suppliers of reusable holloware to hospitals and distributors both nationwide and globally. We offer a comprehensive range manufactured using homo-polymer polypropylene designed for repeated washing and autoclaving at 135C. The "Capecraft" brand has an unrivalled reputation for longevity and durability within the healthcare environment.

Gallipots

Available in ● ●



Pictured L-R: GP80B, GP80R, G60B, GP40B

Lotion Bowls

Available in ● ● ○



Pictured L-R: LB100B (shown with LB100B LID), LB200B, LB150B, LB250B

Kidney Dishes & Vomit Bowl

Available in ● ● ○



Pictured L-R: KD150B, KD200B, KD250B, KD300B & VB12B

Washbowl

Available in ● ○



Pictured: WB350B

Certified to ISO 9001 we are committed to the principals of this internationally recognised quality standard. Our holloware range is designed and manufactured to BS 5452; the British Standard for Hospital Holloware. All of our re-useable holloware can be recycled with polypropylene being identified as a Class 5 material.

Beakers & Feeder Tops

Available in



Pictured: L-R: **BK1N, TFW, TFN**

Funnels

Available in



Pictured: **F100N, F150N**

St. Peter's Boat

Available in



Pictured: **SP1B**

Spatula


Available in



Pictured: **CS1B**
**Please contact sales staff to discuss
Spatula & Bowl Set requirements**

Instrument Tray 300mm

Available in ●●○●●●●●●●●



Pictured: IT300

Instrument Tray 200mm

Available in ●●○●●●●●●●●



Pictured: IT200

Compartmented Instrument Tray 270mm

Available in ●●○●●●●●●●●



Pictured: CIT270

Instrument Tray

Available in ○



Pictured: IT14103N

Instrument Tray Lid

Available in ●



Pictured: IT200B, ITL200B

Instrument Tray with Lid

Available in ●●



Pictured: ITL32B