



Surgical Instrument Trays





Surgical Instrument Trays

Our range of plastic autoclavable sterilisation trays provide a cost-effective system to allow for the decontamination of valuable instruments whilst protecting them during transit; meeting all the specifications required without compromising on quality. Both the base and lid are perforated; the closely fitted lid is transparent to allow for easy identification of the contents. A raised silicone mat is included which has holes to improve drainage during the steam sterilisation process, pips in the tray/mat configuration prevent the mats from sticking. The lightweight trays have integral corner feet to allow for complete circulation during the autoclave process.

- > Lightweight > Plastic > Silicone Mat included
- > Silicone Mats are latex free > Pips in Tray/Mat configuration
- > Perforated:- Tray, Lid & Silicone Mat > Transparent Lid
- > Secure, close fitting Lid > Competitive Price
- > Sterilisation - Steam Autoclaving, Gamma Radiation, Ethylene Oxide
- > Available in 5 sizes/6 configurations > Silicone Inserts available

Reference & Description

	Reference & Description	Length x Width x Depth
ST1A	Micro Tray - Base/Lid/Mat	158mm x 70mm x 19mm
ST2A	Micro Tray - Base/Lid/Mat	203mm x 127mm x 19mm
ST3A	Micro Tray - Base/Lid/Mat	203mm x 127mm x 32mm
ST4A	Micro Tray - Base/Lid/2 Mats/Inner Tray	242mm x 133mm x 38mm
ST4B	Micro Tray - Base/Lid/Mat	242mm x 133mm x 38mm
ST4C	Micro Tray - Base/Lid/Mat	242mm x 133mm x 19mm

